
TR **Vessel Kanül**
Kullanma Kılavuzu

EN **Vessel Cannula**
Instructions For Use

REF: 330 02XX 1

CE 0044

STERILE EO

SCAN FOR INFO

1. Tanım

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EN INSTRUCTIONS FOR USE

1. Description

Bıçakçılar Vessel Cannula is used vascular surgery or cardiopulmonary bypass surgeries to deliver cardioplegia solution or checking the graft integrity. Vessel cannula has PVC and transparent part to insert the vein and a PE material cap to close. The transparent part has been designed with a flexible, tapered tube with atraumatic designed distal tips. It has smooth interior surface to reduce hemolysis. The proximal luer connection is designed to securely connect with a standard cardioplegia delivery line. It is a single use and sterile medical device.

2. Disclaimer of Liability

Bıçakçılar Vessel Cannula is not intended for use other than indicated in the instructions for use. The manufacturer is not responsible for problems arising from inexperience or improper use.

3. Technical Specification

Reference	Check Valve	Length of Body	Cap	Type of Body	Tubing
330 0220 1	without valve	4 cm	Yes	w/o Luer Lock	Flexible
330 0221 1	with valve	4 cm	No	w/o Luer Lock	Flexible
330 0230 1	without valve	5 cm	Yes	w/ Luer Lock	Flexible
330 0231 1	with valve	5 cm	No	w/ Luer Lock	Flexible

4. Mode of Action

Bıçakçılar Vessel Cannula can be used in 2 different ways:

1. By applying an infusion with saline solution, areas leaking from the graft are determined.
2. After completion of the sapheno-coronary anastomosis, infusion of antegrade cardioplegia solutions is performed directly through the graft.

5. Intended Use

Bıçakçılar Vessel Cannula is intended to be used checking graft integrity the vein is gently distended with saline solution to identify leaking segments during harvesting of the "saphenous vein graft". Also, it is intended to be used to facilitate infusion of antegrade cardioplegia solutions directly into the vein graft following completion of the saphenous-coronary anastomosis.

6. Intended User

Bıçakçılar Vessel Cannula is designed to be used by trained healthcare personnel.

7. Intended Environment

Bıçakçılar Vessel Cannula is intended to be used in vascular surgery or cardiopulmonary bypass surgeries.

8. Intended Patient Population

All ages group.

9. Information on Accessories and Other Devices

The Vessel Cannula does not have dedicated accessories. It is used together with compatible medical devices such as cardioplegia delivery systems, syringes, and extracorporeal circulation systems.

10. Material

PVC, PE, Acrylic Polymer, Silicone

11. Indications

Bıçakçılar Vessel Cannula is indicated for use in vascular and cardiac surgical procedures, including coronary artery bypass grafting (CABG), whether performed with or without cardiopulmonary bypass (CPB), to assess vein graft integrity and to facilitate the infusion of cardioplegia or other fluids into the vein graft.

12. Contraindications

Vessel cannulas are not intended for use other than indicated. Vessel cannulas should not be used in patients with known hypersensitivity reactions to the polyvinyl chloride, polyethylene, acrylic polymer, and silicone.

13. Residual Risk

According to the risk analysis carried out by the Bıçakçılar Risk Management Team, the residual risk is determined to be as far as reasonably possible.

14. Recommendations for Use

- Remove the vessel cannula from the packaging using aseptic technique.
- For vessel cannula without valve remove the vented male luer cap and discard it properly.
- After removing the cannula from the packaging, inspect the packaging and product for damage.
- Flush the cannula.
- Using appropriate surgical techniques, introduce the cannula into the vein. Care should be taken during cannula insertion to minimize the possibility of tissue damage by the cannula tip.
- An anchoring suture may be placed proximity to the anchoring hub to secure the vessel onto the cannula.
- Backflush the cannula to remove air.
- Prime the cardioplegia delivery circuit.
- Attach female luer connector of cannula to the cardioplegia set male luer connector.
- After completion of the procedure, remove the cannula using appropriate surgical technique and repair the cannulation site.

Controlling Leakage in a Harvested Vein:

- Connect: Attach the female luer-lock of the vessel cannula to a syringe filled with sterile saline solution.
- Tie Off: Secure one end of the vein with a suture.
- Insert: Place the other end of the vein over one or both different-sized anchor hubs at the distal tip of the vessel cannula.
- Secure: Firmly attach the vein to the cannula.
- Fill: Slowly fill the vein with sterile saline solution.
- Check: Inspect the vein thoroughly for any signs of leakage.

15. Warnings

- Inspect the product for damage prior to use. Damage to the product (i.e., kink crimps, distortions or cuts) may result in obstruction or disruption of flow during use and product performance may be affected. Do not use damaged product.
- If the packaging is damaged, do not use as sterility of the device may have been compromised.
- Use aseptic techniques during the procedure.
- Do not use if protective caps have fallen off or are missing.
- A spare product must always be available during operation.
- Properly affix the vessel cannula to the vein to prevent inadvertent removal during use.
- Do not entrap air during the connection of the vessel cannula to the vein.
- The Lowest possible pressure should be used to infuse cardioplegia solution to minimize hemolysis, air entrapment and, damage to the vessels.
- The product is produced for single use and for single patient use only: the product comes in contact with human blood, body fluids, liquids or gases. Re-sterilization of the product may negatively affect the product performance. Therefore, reuse on other patients may cause cross-contamination, infection and sepsis. Do not resterilize.
- Please monitor the flow during the operation, prevent any complication.

16. Adverse Events and Complications

Complications directly related to the use of vessel cannulas may either be disconnection at the luer connector site or disengagement of the saphenous vein graft from the barbs at the distal end. The cannula shall only be positioned within the intended graft. Improper positioning, including leaving the cannula outside the intended vascular pathway, may result in unintended interaction with the heart and circulatory system, potentially leading to serious complications such as vascular injury, cardiac trauma, and bleeding.

17. Caution

These instructions are provided for informational purposes only. Each surgeon must determine the appropriate use of this device on a case-by-case basis, relying on their medical training, experience, and the specific surgical procedure being performed.

18. Packaging and Transportation

Bıçakcılar Vessel Cannula consists of body, valve, cap assembled and sterile packaged ready to use.

19. Expected Clinical Benefits

Expected clinical benefits of Vessel Cannula are:

- Effective access the vascular system,
- Cannulation the vein graft,
- Reduce leakage complication.

20. Performance Requirements and Limitations

Bıçakcılar Vessel Cannula is used maximum 6 hours during the cardiopulmonary bypass surgeries.

21. Disposal of The Device

The product should be discarded as contaminated medical waste.

22. Notice Regarding Serious Incident

Any serious incident that will occur in relation to the device should be reported to Bıçakcılar and the correspondent competent authority of the Member State in which the user and/or patient is established.

23. Basic UDI-DI

8695476DP330020RP.

24. Shelf Life

The shelf life of the device is 5 years.

25. Storage

The device is recommended to be kept in its original packaging and stored at a temperature between 5°C - 35°C until use. Protect from direct exposure to sunlight. Under proper storage conditions, the device may be used up to the expiry date indicated on the label.

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